

	<b>Citrine</b> "Glistening free"	<b>Aquamarine</b>
A-const (Ultrasound)	119.5	118.0
Optical-const (IOL Master)	119.7	118.3
SRK/T	119.7	118.3
a0 Haigis	1.772	0.846
a1 Haigis	0.400	0.400
a2 Haigis	0.100	0.100
pACD HofferQ	5.97	5.07
SF Holladay 1	2.19	1.29
ACD Holladay 2	5.957	5.140
LF Barrett	2.25	1.52